

Core staff pre-training checklist:

- Assign project to training session in PPMS
 - Change User rights to Novice after the first training for both instruments.
 - Change user rights to autonomous after second training for both instruments.
 - Add lab login to SpectroFlo software (on both instruments) if none exists
 - Request HSRB badge access if needed
-

Two Aurora's (same configuration) available: Winship C-5027 and HSRB E-362B

- Separate badge access request if using the Aurora in Winship (instructions on ECFCC website) <https://sites.winship.emory.edu/facilities>
- First user of the day and instrument is off: warm-up period of 20 minutes is required

Logging in/initializing software

- PPMS
- Launch SpectroFlo software
 - o Login: XXX-Lab
 - o Password: facs123 (
 - o Passwords posted on front of instrument
 - o Acquisition -> opens Acquisition experiment menu

Physical components of cytometer

- Power button
- Visibly check fluid levels on sheath and waste tank (software indicators are not reliable)
 - o Replace sheath when low and/or waste tank is full, leave empty for staff to refill
 - o Spare sheath tanks
 - o Sheath tank – detach line by releasing at the quick disconnect first and then unscrew lid and place on new sheath tank
 - o Waste tank – unscrew, place in holder, pour waste down sink, and add bleach to empty tank

Basics of spectral flow cytometry

- Rather than measuring fluorescence in a single bandpass filter, an entire spectrum of fluorescence is collected for every dye used across a detector array
- 4 laser system (no UV)
- Number of detectors: 48
- Detector arrays off each laser (naming convention – excitation laser, detector #) e.g. YG4

Setting up an Experiment

- Select “New experiment”
- Name Experiment
 - o Experiment name character limit: 50 characters
 - o Include date
- Assigning fluorochromes to experiment
 - o Choose based on laser excitation or type name of fluorochrome to filter
 - o Notify staff to add any fluorochromes not currently in software
- Groups: Creating reference controls for unmixing
 - o Unmixing \neq compensation but achieves the same thing
 - o Autofluorescence extraction
 - Set up an unstained cells reference control for autofluorescence extraction
 - Needs to be the same type of cells in your sample and treated the same way as your sample
 - Group-specific unstained controls
 - o Reference controls are single color spectral controls (same as comp controls)
 - May be cells (recommended) or beads
 - If necessary, add additional unstained controls and assign to each reference control
 - o Labels
- Groups: Experimental samples
 - o Add as many tubes as you have samples
 - o Create cleaning group w/two tubes: rename tubes to bleach and water
- Labels - can be set at experiment, group or tube level
- Keywords - typically not used – for tracking in clinical labs (core staff can add if requested)
- Acquisition
 - o Worksheets (raw or unmixed can be assigned to each group)
 - o Acquisition stopping criteria based on:
 - # of events based on stopping gate
 - Storage gate
 - Stopping time
 - Stopping volume
 - o Number of events to collect for reference controls
 - depends on sample type and marker
 - o 2,000 events minimum in negative and positive populations
 - o For experimental samples, collect enough events of your rarest population to ensure statistical significance
 - o Cleaning tubes
 - Set stopping gate to 1,000,000
 - Set stopping time to 180 seconds (3 minutes)

Acquisition

- Sheath, Waste, Cytometer and Loader status in software
 - o Check cytometer is connected (green)
 - o If using Plate Loader, power should be on (green)
- Raw worksheet orientation
 - o Full spectrum histogram corresponds to individual, color-coded histograms showing individual detectors
 - o Gated vs ungated
- Editing experiment
- Flow rates: low/med/high
 - o uL/min should not be less than 5
 - o Adjusting SIT height/SIT calibration
- Manual SIT flush
- Clean flow cell - stubborn clogs
- Events to display
- Instrument settings
 - o Start w/Cytek Assay settings
 - o SSC detected from two different lasers
 - o Fluorescent detector gains
 - Gains set based on QC every morning
 - Gains are optimized for best resolution (signal to noise ratio)
 - Anything off scale -> titrate Abs and use less dyes
 - Lower fluorescent gains as a last resort
 - Start on excitation laser for fluorochrome being acquired
 - Decrease all detector gains the same amount
 - o Reference controls
 - Always run in Raw worksheet
 - Preview all reference controls first before recording
 - Adjustments to FSC and SSC (more sensitive) and B-SSC will probably be necessary between beads and cells
 - Make any adjustments to fluorescent detector gains
 - SIT extends down into sample during acquisition and retracts otherwise – give it time to retract before removing tube
 - Reference spectral properties on the printout or using Cytek's spectra viewer to confirm the fluorochrome in each reference control

Live unmixing

- Check box for autofluorescence extraction
- Adjust gates for unmixing – adjust in 3 places
 - o Hold “ctrl” down to adjust scatter gate

- Check interval gates includes negatives and positives
- Check correct peak channel is selected (brightest peak on Y axis)
- Select “Live Unmixing” (not Unmix, Save and Open)

Unmixed Worksheets

- Renaming, Save, Save as
- Editing experiment – assigning unmixed worksheet to groups/tubes
- Drawing plots
 - Drawing regions and applying gates
 - Setting acquisition limits based on new gates
 - Creating stats/population hierarchy
 - Editing plot properties – axis scaling
 - Changing parameters on plots
 - AF parameter

Cleaning

- Run on high
- Record tubes
- Water tube stays on SIT
- Check instrument schedule on PPMS
 - Last user of the day performs Fluidics Shutdown
 - 3 mL fluid in each tube
 - Power button on side of instrument

Data

- Two types of FCS files are exported automatically as files are recorded
 - Raw
 - Unmixed
- Shortcut to data folder on Desktop (E drive)
- Edit experiment name in SpectroFlo, not in Windows file explorer
- Use One Drive to transfer data - no flash drives (bookmarked in Chrome)
- Data and experiments are safe for one month

Saving and Templates

- 3 places to save – [add character limits]
 - Experiment – Save as (template) and Save
 - 50 Character Limit for Experiment Name
 - Worksheets
 - 30 Character Limit for user settings Name
 - User Settings – Save as (still be updated based on daily QC)
 - 30 Character Limit for user settings Name

Booking/PPMS

- Emory ID and pw login
- Select instrument, time, then “book”
- Can cancel up to 24 hours in advance of booking for free
- Can move time around on same day
- Requesting assistance – need 48-72 hours advance notice
- Incident reporting
 - o Use low or medium severity
 - o Zoom/remote assistance
 - o Troubleshooting guide

ECFCC website

- Aurora resources
- Cytex webpage (including spectral viewer)