

Laser	Detector	Fluorophores	Initial Voltage Setting	Adj. Voltage Setting	Color	Antibody
Blue (488)	B 515/20	BB515/ A488/ GFP/ FITC				
	B 610/20	BB630				
	B 670/30	BB660				
	B 710/50	PerCP-Cy5.5/ PerCPeF710/ 7-AAD/ PI				
	B 750/30	BB750				
	B 780/60	BB790				
Red (640)	R 670/30	APC/ A633/ A647				
	R 710/50	A700/ APC-R700				
	R 780/60	APC-Cy7/ APC-H7/ Zombie NIR				
Violet (405)	V 431/28	BV421				
	V 470/15	BV510/ Pacific Blue/ BV480				
	V 586/15	BV570/ L/D Aqua/ Pacific Orange				
	V 610/20	BV605				
	V 670/30	BV650				
	V 710/50	BV711				
	V 740/35	BV750				
UV (355)	UV 379/28	BUV395/ Indo-1 Bound				
	UV 515/30	BUV496/ Hoechst 33342/ Indo-1 Unbound				
	UV 586/15	BUV 563				
	UV 670/30	BUV661				
	UV 740/35	BUV737				
	UV 820/60	BUV805				
	Yellow/ Green (561)	YG 586/15	PE/ BYG584			
YG 610/20		PE-CF594/ RFP				
YG 670/30		PE-Cy5/ 7-AAD				
YG 710/50		PE-Cy5/ 7-AAD				
	YG 780/60	PE-Cy7/ BYG790				